


General Bed
25000 EB100







Technical Specifications

Overall dimensions	2310 mm × 1080 mm
Mattress platform	2070 mm × 800 mm
Height	560 mm
Safety side height	390 mm
Back rest angle	0 - 85°
Knee rest angle	0 - 35°
Bed weight	102 kg
Safe working load	250 kg

Casters	125 mm
Power input max.	250 VA
Main power system	220 V – 50 Hz, 24V DC
Electrical current protection	Class I
Electrical shock protection	Type B
Operation mode	10%, max.: 2 min/18 min
IP rating	IP 66

Technical and Quality Assurance Standards
 IEC 60601- 2 - 52
 ISO 13485
 ISO 9001
 ISO 14971

Positions



Features



Safety Sides

25000 EB100 has blow molded retractable split safety sides that are equipped with embedded control keyboards for the patient and caregiver. Providing a protective height and length, these parallel safety sides adhere to the highest standards of safety to prevent patient fall and entrapment.



Embedded Control Panel

User friendly embedded control panels are located on both safety sides of head section. Using the control keyboards on internal and external parts of safety sides, the patient and caregiver are able to adjust the back rest and knee rest positions of the bed.



Bumpers

Rotary bumpers that are located on four corners of the bed can prevent to a large extent the damages that may occur due to unwanted collision of the bed to the surroundings during transportation.



Urine Bag Holder

Urine bag holders can be easily installed on either side of the bed in one step.



Bed Ends

The seamless blow molded bed ends of 25000 EB100 are lightweight, durable, and easy to clean. The bed ends are of high quality and compatible for steam cleaning. In order to have a better access to patient's head in emergency cases, the caregiver can easily remove the bed ends by grasping the hand holds and lifting them up in one single step. Moreover, the bed ends have integrated steel fittings which make them extremely rigid; so the nursing staff can move the bed by holding and pushing them without worries.




Casters

125 mm casters run extremely smoothly; the quality that makes it significantly easy for the nursing staff to transport the bed in the medical centers. Having diagonal locks, the casters can be locked in order to prevent the bed from unintended movement.



Mattress Platforms

The ABS mattress platforms of this general bed are produced by injection molding method. These platforms have some remarkable features such as smooth surfaces, high durability, and easy to remove for cleaning and disinfection. Plus, the mattress platforms are compatible for steam cleaning.



Actuators

All positions of this bed are achieved by powerful, reliable, and low noise actuators. As the actuators have IP66 rating, they are protected against ingress of dust and water. This prominent feature of actuators makes it easy to clean and disinfect the bed in the medical centers on a regular basis.