

General Bed Model 25000MD1



Technical Data

Overall dimensions	2255 mm L x 1048 mm W
Mattress platform	2000 mm L x 810 mm W
Height	580 mm
Side rail level above mattress platform	415 mm
Back rest angle	67°
Knee rest angle	34°
Bed weight	95 kg
Safe working load	220 kg
Castors (With two independent brakes)	125 mm

Standards

IEC 60601-2-52
ISO 13485
ISO 9001
ISO 14971

Features

- Back rest & knee rest positions can be adjusted by two cranks
- Removable blow molded head & foot boards with steel supports;
- one holder in each corner to installing accessories (IV stand, Traction bar, etc.)
- A unique and stylish design of parallel rail board combines blow molded plastic parts and steel bars to provide a protective height and length. The rail board adhere to the highest standards of safety to prevent patient fall and entrapment.
- Four removable injected ABS mattress platforms
- 125 mm diameter high quality castors with two independent brakes
- Detachable urine bag holders on both sides
- Capability of installing accessories (Orthopedic traction frame)
- Phosphated pretreatment for coating

