

General Bed Model 25000EC1S



Technical Data

■ Overall dimensions	2250 mm L x 1010 mm W
■ Mattress platform	2000 mm L x 810 mm W
■ Height	540 mm
■ Side rail level above mattress platform	420 mm
■ Back rest angle	78°
■ Knee rest angle	30°
■ Bed weight	87 kg
■ Safe working load	250 kg
■ Castors (With two independent brakes)	150 mm

Standards

IEC 60601-2-52
ISO 13485
ISO 9001
ISO 14971

Features

- Back rest & knee rest positions can be adjusted by two actuators
- Removable blow molded head & foot boards with steel supports;
- one holder in each corner to installing accessories (IV stand, Traction bar, etc.)
- A unique and stylish design of parallel rail board combines blow molded plastic parts and steel bars to provide a protective height and length. The rail board adhere to the highest standards of safety to prevent patient fall and entrapment.
- Four removable injected ABS mattress platforms
- 150 mm diameter high quality castors with two independent brakes
- Detachable urine bag holders on both sides
- Capability of installing accessories (Orthopedic traction frame, Lifting pole)
- Phosphated pretreatment for coating
- Electrostatic powder coating

